

# Small Animal I.V. Set One Time Purchases

## Shoe Mounting Kit

Reorder No. SH9100

The shoe and tubing tabs hold the swivel device and primary I.V. set tubing to the top of the cage. The shoe and tubing tabs have adhesive backing for permanent attachment to the top of the cage. The tubing and swivel device are snapped into place during in-cage I.V. treatment.



Kit includes alcohol prep pads for cleaning the top of the cage before affixing shoe and tubing tabs.

Swivel Device Mounting Shoe




Tubing Tabs



## Arm Strap With Shoe

Reorder No. SH9200

Blue carrying strap designed for comfort and convenience. It permits one person to carry or walk a patient, along with the fluids, the Flex-coil® Swivel set, and I.V. tubing, while the patient continues to receive I.V. therapy. Shoe has an attached "D" ring with a Velcro-backed strap that can be detached from the arm strap to temporarily hang the fluids, the Flex-coil® Swivel device and tubing in another area. The arm strap is a must to improve hospital efficiency.



Velcro-backed Detachable Strap

Fluid Bag Tabs

Swivel Device Mounting Shoe

"D" Ring



## INTERNATIONAL WIN, LIMITED

340 North Mill Road, Suite 6, Kennett Square, PA 19348-2853 USA • (610) 444-0170 • (800) 359-4946 • Fax (610) 444-0171  
e-mail: [info@internationalwin.com](mailto:info@internationalwin.com) • [www.internationalwin.com](http://www.internationalwin.com)