

Set-up Instructions - Flex-coil® Swivel System on Cats and Dogs

U.S. Patent No. 5,163,923 • Canadian Patent No. 2,119,376

STEP FIVE

Prime the appropriate sized Extension Set and close the on/off clamp. The 14 inch Extension Set with KINK-SAFE™ Adapter (EXT14) is ideal for cats, small dogs and puppies; the 24 inch Extension Set with KINKSAFE™ Adapter (EXT24) is ideal for average sized dogs; and, the 36 inch Extension Set with KINK-SAFE™ Adapter (EXT36) is ideal for giant breed dogs. Attach the male Luer connector of the KINK-SAFE™ adapter to the hub of the I.V. catheter. With a syringe attached to the female Luer connector on the Extension Set flush several milliliters of fluid through the Extension Set and I.V. catheter to check the patency of the vein.



STEP SIX

Attach the male Luer connector on the end of the coil to the female Luer connector at the proximal end of the Extension Set.

Depress the tab on the releasable cable tie at the end of the coil, open the tie, pass it around a collar, releasable collar, identification band or harness and secure the cable tie so that end of the coil is positioned behind the head between the shoulder blades (out of sight and out of easy reach).



STEP SEVEN

The final steps are to open the on/off clamp on the Extension Set and set the rate of fluid flow on the infusion pump or with the flow control device on the primary I.V. set.



INTERNATIONAL WIN, LIMITED

340 North Mill Road, Suite 6, Kennett Square, PA 19348-2853 USA • (610) 444-0170 • (800) 359-4946 • Fax (610) 444-0171
e-mail: info@internationalwin.com • www.internationalwin.com