

**FOR IMMEDIATE RELEASE**  
**March 17, 2009**

**I**NTERNATIONAL WIN, LTD., Kennett Square, PA has completed its testing and recertification of the revised and improved Flex-coil® Swivel System line of products. The Flex-coil® Swivel System line of fluid I.V. Sets has been rigorously tested and has been found to be fully functional and approved for sale commencing on February 10, 2009. The Flex-coil® Swivel System line of products includes the CR9003 three foot, CR9005 five foot and CR9010 ten foot coiled Core Flex-coil® Swivel fluid delivery sets with internal swivels; and, the PCRE9005, Standard Flex-coil Kit and the PMCRE9005, Pediatric Flex-coil® Swivel Kit with five foot coils and internal swivels.

The basic purpose, design and function of the Flex-coil® Swivel System line of products have not changed. The improvements to the Flex-coil® Swivel Sets include a new K-seal to better protect the devices from leaks of air and water during use, application of a silicone lubricant to reduce the seal to pin friction which improves ease of pin rotation, a redesigned central rotating pin with tubing barb, and addition of a crimp ring to more securely attach the polyurethane coil to the rotating central pin.

During production of the Flex-coil® Swivel devices, additional steps have been introduced to reduce the likelihood of air and water leaks during use and to insure the central pin rotates freely when first placed in service. An elaborate and sensitive test system has been designed to check every swivel device for minute air and water leaks and if leaks are revealed the swivel is discarded. Each central pin is also rotated at high speed for 30 seconds to insure the silicone lubricant is distributed evenly and completely between the seal and central rotating pin.

For further information contact:

William J. Donawick, D.V.M.  
International WIN, Ltd.  
340 North Mill Road  
Suite 6  
Kennett Square, PA 19348-2853 U.S.A.  
610-444-0170  
800-359-4946

