


Small Animal Swivel Run Rod

Reorder No. SRR9000

Adjustable length, portable, 1 inch diameter tension rod with attached Shoe Mount and tubing clips for use in runs. Rod with shoe suspends the Flex-coil® Swivel Set above and out of reach of the patient at the center of the run and enables the patient to move freely without kinking, twisting or damaging the coil. The rod can be mounted on steel or masonry walls, chain-link fencing, and horizontal or vertical bar side walls. Swivel Run Rods are made to order to best fit your runs (call us for details).



INTERNATIONAL WIN, LIMITED

340 North Mill Road, Suite 6, Kennett Square, PA 19348-2853 USA • (610) 444-0170 • (800) 359-4946 • Fax (610) 444-0171
e-mail: info@internationalwin.com • www.internationalwin.com